

Anti-Cytokeratin 6 (Krt6) Rabbit Polyclonal Antibody For Wb And Elisa - P02538

Item Number: CM0030483



Introduction

Rabbit polyclonal antibody specific to human Cytokeratin 6 (KRT6), validated for WB and ELISA. Detects a 60 kDa band; ideal for studying epithelial integrity, wound repair, and keratinopathies. - P02538

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 6 (KRT6) Rabbit pAb
Remarks/Alias	K6A; K6C; K6D; PC3; CK6A; CK6C; CK6D; CK-6C; CK-6E; KRT6C; KRT6D; Cytokeratin 6 (KRT6)
Species	Human
Gene ID (Human)	3853
Gene ID	3853
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cytokeratin 6 (KRT6) (NP_005545.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02538
Protein Weight	60kDa
Shipping	Ice bag