

Anti-Cytokeratin 7 (Krt7) Monoclonal Antibody For Wb, If, Ihc, Elisa - P08729

Item Number: CM0005232



Introduction

High-quality rabbit monoclonal antibody targeting human Cytokeratin 7 (KRT7, CK7). Appropriate for WB, IF, IHC, ELISA detection. Cross-reacts with mouse and rat orthologs. Recombinant immunogen (aa 291-463, P08729).

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 7 (KRT7) Rabbit mAb
Remarks/Alias	K7; KRT7; SCL; K2C7; Cytokeratin 7 (CK7)
Species	Human
Gene ID (Human)	3855
Gene ID	3855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 291-463 of human Cytokeratin 7 (KRT7) (P08729)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08729
Protein Weight	51 kDa
Shipping	Ice bag