

Anti-Cytokeratin 7 (Krt7) Monoclonal Antibody For Wb, Ihc-P, If, Elisa - P08729

Item Number: CM0004568



Introduction

Rabbit monoclonal anti-Cytokeratin 7 (KRT7) antibody validated for WB, IHC-P, IF, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for epithelial biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 7 (KRT7) Rabbit mAb
Remarks/Alias	K7; CK7; SCL; K2C7; Cytokeratin 7 (KRT7)
Species	Human
Gene ID (Human)	3855
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cytokeratin 7 (KRT7) (P08729)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08729
Protein Weight	51kDa
Shipping	Ice bag