

Anti-Cytokeratin 7 (Krt7) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P08729

Item Number: CM0024480



Introduction

Cytokeratin 7 (KRT7) rabbit polyclonal antibody validated for WB, IF-P, ELISA. Reacts with human, mouse, rat. Targets a cytoskeletal protein involved in interferon signaling and viral translational regulation.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 7 (KRT7) Rabbit pAb
Remarks/Alias	K7; CK7; SCL; K2C7; Cytokeratin 7 (KRT7)
Species	Human
GeneID (Human)	3855
GeneID	3855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-469 of human Cytokeratin 7 (KRT7) (NP_005547.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08729
Protein Weight	51kDa
Shipping	Ice bag