

## Anti-Cytokeratin 8 (Ck8) Monoclonal Antibody For Wb And Elisa - P05787

Item Number: CM0004802



### Introduction

Mouse monoclonal antibody targeting human Cytokeratin 8 (CK8) for WB and ELISA applications. Validated for human, mouse, and rat species, with epitope mapping to a synthetic peptide. Ideal for epithelial cell and keratin research.

[Learn More](#)

| Parameter        | Value                                               |
|------------------|-----------------------------------------------------|
| Product Name     | Cytokeratin 8 Mouse mAb                             |
| Remarks/Alias    | K8; KO; CK8; CK-8; CYK8; K2C8; CARD2; Cytokeratin 8 |
| Species          | Human                                               |
| GeneID           | 3856                                                |
| Immunogen        | Synthetic peptide of CK8                            |
| Source           | Mouse                                               |
| Category         | Monoclonal Antibodies                               |
| Application      | WB, ELISA                                           |
| Cross Reactivity | Mouse, Rat                                          |
| SWISS            | P05787                                              |
| Protein Weight   | 54 kDa/56 kDa                                       |
| Shipping         | Ice bag                                             |