


Anti-Cytokeratin 8 (Krt8) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05787

Item Number: CM0016189



Introduction

Rabbit polyclonal antibody targeting Cytokeratin 8 (KRT8). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for studying keratin intermediate filaments and epithelial cell mechanics.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 8 (KRT8) Rabbit pAb
Remarks/Alias	K8; KO; CK8; CK-8; CYK8; K2C8; CARD2; Cytokeratin 8 (KRT8)
Species	Human
Gene ID (Human)	3856
Gene ID	3856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-483 of human Cytokeratin 8 (KRT8) (NP_002264.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05787
Protein Weight	54 kDa/56 kDa
Shipping	Ice bag