

Anti-Cytokeratin 9 (Krt9) Rabbit Monoclonal Antibody For Wb, Elisa - P35527

Item Number: CM0010148

Introduction

Rabbit monoclonal antibody targeting human Cytokeratin 9 (KRT9), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for studying palmoplantar epidermis differentiation and related keratin disorders.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 9 (KRT9) Rabbit mAb
Remarks/Alias	K9; CK-9; EPPK; Cytokeratin 9 (KRT9)
Species	Human
GeneID (Human)	3857
GeneID	3857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 254-412 of human Cytokeratin 9 (KRT9) (P35527)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35527
Protein Weight	62kDa
Shipping	Ice bag