

Anti-Cytokeratin 9 (Krt9) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P35527

Item Number: CM0016100

Introduction

High-quality rabbit polyclonal antibody against human Cytokeratin 9 (KRT9), validated for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat KRT9. Ideal for palmoplantar keratinocyte differentiation research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 9 (KRT9) Rabbit pAb
Remarks/Alias	K9; CK-9; EPPK; Cytokeratin 9 (KRT9)
Species	Human
GeneID (Human)	3857
GeneID	3857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-450 of human Cytokeratin 9 (KRT9) (NP_000217.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35527
Protein Weight	62kDa
Shipping	Ice bag