


Anti-Cytokeratin Ae1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P13645 / P02533 / P19012 / P08779 / P08727

Item Number: CM0029508



Introduction

Cytokeratin AE1 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Recognizes human, mouse, rat keratins (KRT10/14/15/16/19) with molecular weights 44-59 kDa. Ideal for epithelial differentiation and carcinoma research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin AE1 Rabbit pAb
Remarks/Alias	BIE; EHK; IHL; K10; KPP; BCIE; CK10
Species	Human
GeneID (Human)	3858/3861/3866/3868/3880
GeneID	3858/3861/3866/3868/3880
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13645/P02533/P19012/P08779/P08727
Protein Weight	59kDa/52kDa/49kDa/51kDa/44kDa
Shipping	Ice bag