


Anti-Cytokeratin Ae1 Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P13645 / P02533 / P19012 / P08779 / P08727

Item Number: CM0008515



Introduction

Rabbit monoclonal antibody against Cytokeratin AE1, detecting multiple type I keratins (K10, K14, K15, K16, K19). Validated for IF-P, IHC-P, and ELISA in human, mouse, and rat tissues. Ideal for epithelial cell biology and carcinoma research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin AE1 Rabbit mAb
Remarks/Alias	BIE; EHK; IHL; K10; KPP; BCIE; CK10
Species	Human
GeneID (Human)	3858/3861/3866/3868/3880
GeneID	3858/3861/3866/3868/3880
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13645/P02533/P19012/P08779/P08727
Protein Weight	59kDa/52kDa/49kDa/51kDa/44kDa
Shipping	Ice bag