


Anti-Cytokeratin Ae3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04264/P35908/P12035/P19013/P13647/P02538/P08729

Item Number: CM0025134



Introduction

Rabbit polyclonal antibody targeting Cytokeratin AE3 (K1, CK1, EHK) validated for Western blot, IHC-P, and ELISA. Reacts with human samples. Ideal for epithelial and keratinocyte biology research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin AE3 Rabbit pAb
Remarks/Alias	K1; CK1; EHK; AEI2; EHK1; EPPK; KRT1A; NEPPK
Species	Human
GeneID (Human)	3848/3849/3850/3851/3852/3853/3855/3856
GeneID	3848/3849/3850/3851/3852/3853/3855/3856
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 181-280 of human Cytokeratin AE3 (NP_006112.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P04264/P35908/P12035/P19013/P13647/P02538/P08729
Protein Weight	66kDa/65kDa/64kDa/56kDa/62kDa/60kDa/51kDa
Shipping	Ice bag