


Anti- δ -Catenin/p120 Catenin Polyclonal Antibody For Wb, Ihc-P, Elisa - O60716

Item Number: CM0023865



Introduction

Rabbit polyclonal antibody targeting δ -Catenin/p120 Catenin (CTNND1) for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat. Key regulator of cadherin-mediated adhesion and transcription; ideal for cancer and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	δ -Catenin/p120 Catenin Rabbit pAb
Remarks/Alias	CAS; p120; BCDS2; CTNND; P120CAS; P120CTN; p120(CAS); p120(CTN); δ -Catenin/p120 Catenin
Species	Human
GeneID (Human)	1500
GeneID	1500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 573-832 of human δ -Catenin/p120 Catenin (NP_001078938.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60716
Protein Weight	108kDa
Shipping	Ice bag