

Anti-D2Hgdh Polyclonal Antibody For Wb And Elisa - Q8N465

Item Number: CM0024032



Introduction

Rabbit polyclonal antibody targeting human D2HGDH for WB and ELISA applications. Cross-reacts with mouse and rat. Ideal for studying D-2-hydroxyglutarate metabolism, mitochondrial function, and related disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	D2HGDH Rabbit pAb
Remarks/Alias	D2HGD; D2HGDH
Species	Human
Gene ID	728294
Immunogen	Recombinant fusion protein containing amino acids 322-521 of human D2HGDH (NP_689996.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N465
Protein Weight	56kDa
Shipping	Ice bag