


Anti-D53 Polyclonal Antibody For Wb And Elisa - Q2Rbp2

Item Number: CM0022897



Introduction

Rabbit polyclonal antibody targeting DWARF 53 (D53) from *Oryza sativa*. Validated for WB and ELISA. Useful for studying strigolactone signaling and tiller development in rice. Molecular weight approximately 123 kDa.

[Learn More](#)

Parameter	Details
Product Name	D53 Rabbit pAb
Remarks/Alias	D53; DWARF53; Osj_32651
Species	<i>Oryza sativa</i>
GeneID	4349543
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1000-1100 of <i>Oryza sativa</i> D53 (Q2RBP2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q2RBP2
Protein Weight	123kDa
Shipping	Ice bag