

Anti-Daam2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q86T65

Item Number: CM0025387

Introduction

Rabbit polyclonal antibody targeting DAAM2, suitable for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Essential tool for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	DAAM2 Rabbit pAb
Remarks/Alias	NPHS24; dj90A20A.1; DAAM2
Species	Human
GeneID (Human)	23500
GeneID	23500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human DAAM2 (NP_056160.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86T65
Protein Weight	123kDa
Shipping	Ice bag