

Anti-Dab1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75553

Item Number: CM0016299

Introduction

Rabbit polyclonal DAB1 antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Targeting the reelin signaling adapter, ideal for neuronal migration and brain development studies.

[Learn More](#)

Parameter	Value
Product Name	DAB1 Rabbit pAb
Remarks/Alias	SCA37; DAB1
Species	Human
Gene ID (Human)	1600
Gene ID	1600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human DAB1 (NP_066566.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75553
Protein Weight	64kDa
Shipping	Ice bag