


Anti-Dab2 Rabbit Polyclonal Antibody For Wb, Fc, Elisa - P98082

Item Number: CM0028682



Introduction

High-quality rabbit polyclonal antibody targeting human disabled homolog 2 (DAB2/DOC2), validated for WB, FC, ELISA. Key reagent for endocytosis, signaling, and cancer research. Specific to human DAB2.

[Learn More](#)

Parameter	Value
Product Name	DAB2 Rabbit pAb
Remarks/Alias	DOC2; DOC-2; DAB2
Species	Human
Gene ID (Human)	1601
Gene ID	1601
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 340-580 of human DAB2 (NP_001334.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P98082
Protein Weight	82kDa
Shipping	Ice bag