


Anti-Dach1 Polyclonal Antibody For Wb, Elisa - Q9Ui36

Item Number: CM0024781



Introduction

Rabbit polyclonal antibody against human DACH1, validated for WB and ELISA. Reacts with human, mouse, rat samples. Detects 79 kDa protein, a transcription factor involved in organogenesis and TGF-beta signaling.

[Learn More](#)

Parameter	Value
Product Name	DACH1 Rabbit pAb
Remarks/Alias	DACH; DACH1
Species	Human
Gene ID (Human)	1602
Gene ID	1602
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Ui36
Protein Weight	79kDa
Shipping	Ice bag