

Anti-Dact3 Rabbit Polyclonal Antibody For Wb And Elisa - Q96B18

Item Number: CM0030035

Introduction

CamelBio's DACT3 rabbit polyclonal antibody targets the N-terminus (aa1-160) of human DACT3, recognizing mouse and rat DACT3. Validated for WB and ELISA, it is ideal for Wnt signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	DACT3 Rabbit pAb
Remarks/Alias	RRR1; DAPPER3; DACT3
Species	Human
GenelD (human)	147906
GenelD	147906
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96B18
Protein Weight	65kDa
Shipping	Ice bag