

Anti-Dad1 Rabbit Polyclonal Antibody (Wb, Elisa) - P61803

Item Number: CM0018496



Introduction

DAD1 Rabbit Polyclonal Antibody validated for WB, ELISA. Detects human, mouse, rat DAD1, subunit of the OST complex crucial for N-glycosylation. Ideal for ER stress and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	DAD1 Rabbit pAb
Remarks/Alias	OST2; DAD1
Species	Human
GeneID (Human)	1603
GeneID	1603
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DAD1 (NP_001335.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61803
Protein Weight	12kDa
Shipping	Ice bag