

Anti-Dag1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q14118

Item Number: CM0024029



Introduction

Rabbit polyclonal antibody against human DAG1, validated for WB, IF/ICC, and ELISA. Detects dystroglycan in human, mouse, and rat. Ideal for studying cell adhesion, muscular dystrophy, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	DAG1 Rabbit pAb
Remarks/Alias	A3a; DAG; AGRNR; 156DAG; MDDGA9; MDDGC7; MDDGC9; LGMDR16; DAG1
Species	Human
GeneID (human)	1605
GeneID	1605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 654-749 of human DAG1 (NP_004384.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14118
Protein Weight	97kDa
Shipping	Ice bag