

Anti-Dag1 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q14118

Item Number: CM0016072



Introduction

Rabbit polyclonal antibody against human DAG1 (Dystroglycan 1), validated for WB, IF/ICC, IP, ELISA. Suitable for dystroglycan complex signaling, extracellular matrix and receptor studies.

[Learn More](#)

Parameter	Value
Product Name	DAG1 Rabbit pAb
Remarks/Alias	A3a; DAG; AGRNR; 156DAG; MDDGA9; MDDGC7; MDDGC9; LGMDR16; DAG1
Species	Human
Gene ID	1605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 776-895 of human DAG1 (NP_001171108.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14118
Protein Weight	97kDa
Shipping	Ice bag