

Anti-Dao Monoclonal Antibody For Wb, Elisa - P14920

Item Number: CM0005933



Introduction

Rabbit monoclonal antibody against human D-amino-acid oxidase (DAO), validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for studying D-amino acid metabolism and neurological pathways.

[Learn More](#)

Parameter	Value
Product Name	DAO Rabbit mAb
Remarks/Alias	DAAO; OXDA; DAMOX; DAO
Species	Human
GeneID (Human)	1610
GeneID	1610
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DAO (NP_001908.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14920
Protein Weight	39kDa
Shipping	Ice bag