

Anti-Dao Rabbit Polyclonal Antibody For Wb, If, Elisa - P14920

Item Number: CM0025822



Introduction

Rabbit polyclonal antibody against human D-amino-acid oxidase (DAO). Validated for WB, IF, and ELISA applications in human, mouse, and rat samples. Ideal for metabolism, neuroscience, and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	DAO Rabbit pAb
Remarks/Alias	DAAO; OXDA; DAMOX; DAO
Species	Human
Gene ID (Human)	1610
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-347 of human DAO (NP_001908.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14920
Protein Weight	39kDa
Shipping	Ice bag