

Anti-Dap Kinase 1 (Dapk1) Monoclonal Antibody For Wb And Elisa - P53355

Item Number: CM0002464



Introduction

Rabbit monoclonal antibody targeting DAP Kinase 1 (DAPK1) residues 250-350. Suitable for WB and ELISA; cross-reacts with human and mouse. For research in apoptosis, autophagy, and cell survival signaling.

[Learn More](#)

Parameter	Value
Product Name	DAP Kinase 1 (DAPK1) Rabbit mAb
Remarks/Alias	DAPK, ROCO3, DAP Kinase 1 (DAPK1)
Species	Human
Gene ID	1612
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P53355
Protein Weight	160kDa
Shipping	Ice bag