


Anti-Dap3 Rabbit Polyclonal Antibody For Wb, Elisa - P51398

Item Number: CM0025719



Introduction

Rabbit polyclonal antibody against human DAP3 (MRPS29), validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for mitochondrial biology and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	DAP3 Rabbit pAb
Remarks/Alias	DAP-3; S29mt; MRPS29; MRP-S29; bMRP-10; DAP3
Species	Human
Gene ID (Human)	7818
Gene ID	7818
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 159-398 of human DAP3 (NP_387506.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51398
Protein Weight	46kDa
Shipping	Ice bag