


Anti-Dapk3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43293

Item Number: CM0013285



Introduction

Rabbit polyclonal antibody targeting DAPK3 (ZIPK/DLK). Suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Supports research on apoptosis, autophagy, and actin dynamics.

[Learn More](#)

Parameter	Value
Product Name	DAPK3 Rabbit pAb
Remarks/Alias	DLK; ZIP; ZIPK; DAPK3
Species	Human
GeneID (human)	1613
GeneID	1613
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DAPK3 (NP_001339.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43293
Protein Weight	53kDa
Shipping	Ice bag