

Anti-Dapp1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Un19

Item Number: CM0013035



Introduction

Rabbit polyclonal antibody against human DAPP1 (BAM32), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for B-cell receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	DAPP1 Rabbit pAb
Remarks/Alias	BAM32; DAPP1
Species	Human
GeneID (Human)	27071
GeneID	27071
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 100-273 of human DAPP1 (NP_055210.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UN19
Protein Weight	32kDa
Shipping	Ice bag