

Anti-Darpp32 Monoclonal Antibody For Wb, If/Icc, If-F, Elisa - Q9Ud71

Item Number: CM0007325

Introduction

Validated rabbit monoclonal anti-DARPP32 antibody for WB, IF/ICC, IF-F, and ELISA. Detects DARPP32 in human, mouse, and rat, ideal for dopamine signaling and PP1 inhibition studies.

[Learn More](#)

Parameter	Value
Product Name	DARPP32 Rabbit mAb
Remarks/Alias	DARPP32; DARPP-32
Species	Human
Gene ID (Human)	84152
Gene ID	84152
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UD71
Protein Weight	23 kDa
Shipping	Ice bag