

Anti-Darpp32 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ud71

Item Number: CM0008237

Introduction

Recombinant rabbit monoclonal antibody targeting human DARPP32 (PPP1R1B, DARPP-32). Validated for WB and ELISA, cross-reacts with mouse and rat. A key tool for studying dopamine/cAMP-regulated protein phosphatase 1 inhibition and neuronal signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	DARPP32 Rabbit mAb
Remarks/Alias	DARPP32; DARPP-32
Species	Human
GeneID (human)	84152
GeneID	84152
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UD71
Protein Weight	23kDa
Shipping	Ice bag