

Anti-Darpp32 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ud71

Item Number: CM0024652

Introduction

Rabbit polyclonal antibody against human DARPP32 (PPP1R1B), a dopamine- and cAMP-regulated neuronal phosphoprotein inhibitor of PP1. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Predicted band ~23 kDa.

[Learn More](#)

Parameter	Value
Product Name	DARPP32 Rabbit pAb
Remarks/Alias	DARPP32; DARPP-32
Species	Human
Gene ID (Human)	84152
Gene ID	84152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-116 of human DARPP32 (NP_115568.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UD71
Protein Weight	23kDa
Shipping	Ice bag