

Anti-Dars Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P14868

Item Number: CM0029330



Introduction

Rabbit polyclonal antibody targeting human DARS (Aspartate-tRNA ligase). Suitable for WB, IHC-P, IP, and ELISA. Detects endogenous levels in human, mouse, and rat samples. Ideal for studies on protein translation, amino acid metabolism, and neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	DARS Rabbit pAb
Remarks/Alias	DARS; HBSL; aspRS
Species	Human
Gene ID (Human)	1615
Gene ID	1615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 222-501 of human DARS (NP_001340.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14868
Protein Weight	57kDa
Shipping	Ice bag