

Anti-Dars2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Pi48

Item Number: CM0030223



Introduction

DARS2 rabbit polyclonal antibody targeting mitochondrial aspartyl-tRNA synthetase. Suitable for Western blot, IHC-P, and ELISA in human and mouse samples. Useful for mitochondrial translation research.

[Learn More](#)

Parameter	Value
Product Name	DARS2 Rabbit pAb
Remarks/Alias	LBSL; ASPRS; mtAsPRS; MT-ASPRS; DARS2
Species	Human
GeneID (Human)	55157
GeneID	55157
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6PI48
Protein Weight	74kDa
Shipping	Ice bag