

Anti-Daxx Rabbit Monoclonal Antibody For Wb, Elisa - Q9Uer7

Item Number: CM0003119



Introduction

Anti-Daxx rabbit monoclonal antibody validated for human samples in WB and ELISA. Targets death domain-associated protein 6, a transcription corepressor and histone chaperone involved in apoptosis and chromatin remodeling. SWISS: Q9UER7.

[Learn More](#)

Parameter	Value
Product Name	Daxx Rabbit mAb
Remarks/Alias	DAP6; EAP1; BING2; SMIM40; Daxx
Species	Human
GeneID (Human)	1616
GeneID	1616
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Daxx (Q9UER7))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UER7
Protein Weight	81kDa
Shipping	Ice bag