

Anti-Daxx Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uer7

Item Number: CM0023141



Introduction

This Daxx rabbit polyclonal antibody detects human and mouse Daxx, a death domain-associated protein 6. Validated for WB and ELISA, it is suitable for research on transcriptional regulation, apoptosis, chromatin remodeling, and viral infection.

[Learn More](#)

Parameter	Value
Product Name	Daxx Rabbit pAb
Remarks/Alias	DAP6; EAP1; BING2; SMIM40; Daxx
Species	Human
GeneID (Human)	1616
GeneID	1616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-400 of human Daxx (NP_001135441.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UER7
Protein Weight	81kDa
Shipping	Ice bag