

Anti-Dazap1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96Ep5

Item Number: CM0003127



Introduction

DAZAP1 Rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, and rat DAZAP1. Ideal for reproductive biology and spermatogenesis research.

[Learn More](#)

Parameter	Value
Product Name	DAZAP1 Rabbit mAb
Species	Human
Gene ID (Human)	26528
Immunogen	Recombinant fusion protein containing amino acids 265-407 of human DAZAP1 (Q96EP5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EP5
Protein Weight	43kDa
Shipping	Ice bag