

## Anti-Dazap2 Rabbit Polyclonal Antibody For Wb And Elisa - Q15038

Item Number: CM0017925

### Introduction

Rabbit polyclonal antibody against human DAZAP2, validated for WB and ELISA. Cross-reacts with mouse. Ideal for studying HIPK2 degradation, DNA damage response, and Wnt signaling.

[Learn More](#)

| Parameter                  | Value                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name               | DAZAP2 Rabbit pAb                                                                                                 |
| Remarks/Alias              | PRTB; DAZAP2                                                                                                      |
| Species                    | Human                                                                                                             |
| Gene ID (Human)            | 9802                                                                                                              |
| Gene ID                    | 9802                                                                                                              |
| Immunogen                  | Recombinant fusion protein containing a sequence corresponding to amino acids 1-168 of human DAZAP2 (NP_055579.1) |
| Source                     | Rabbit                                                                                                            |
| Category                   | Polyclonal Antibodies                                                                                             |
| Application                | WB, ELISA                                                                                                         |
| Cross Reactivity           | Mouse                                                                                                             |
| UniProt Accession          | Q15038                                                                                                            |
| Predicted Molecular Weight | 17kDa                                                                                                             |
| Shipping Condition         | Ice bag                                                                                                           |