

Anti-Dazl Rabbit Polyclonal Antibody For Wb, Elisa - Q92904

Item Number: CM0012982

Introduction

Rabbit polyclonal antibody against human DAZL, validated for WB and ELISA. Cross-reacts with human and mouse. DAZL is an RNA-binding protein essential for gametogenesis, making this antibody ideal for reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	DAZL Rabbit pAb
Remarks/Alias	DAZH; DAZL1; DAZLA; SPGYLA; DAZL
Species	Human
GeneID (Human)	1618
GeneID	1618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 186-295 of human DAZL (NP_001342.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92904
Protein Weight	33kDa
Shipping	Ice bag