


Anti-Dbh Polyclonal Antibody For Wb, Ihc-P, Elisa - P09172

Item Number: CM0018626



Introduction

Rabbit anti-DBH polyclonal antibody validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat dopamine beta-hydroxylase. Ideal for neurotransmitter biosynthesis studies. Ships on ice bag.

[Learn More](#)

Parameter	Value
Product Name	DBH Rabbit pAb
Remarks/Alias	DBM, ORTHYP1, DBH
Species	Human
GeneID (Human)	1621
GeneID	1621
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 388-617 of human DBH (NP_000778.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09172
Protein Weight	69kDa
Shipping	Ice bag