


Anti-Dbh Rabbit Monoclonal Antibody For Wb, Ip, If, Ihc-P, Elisa - P09172

Item Number: CM0010758



Introduction

Rabbit monoclonal antibody against human DBH, validated for WB, IP, IF, IHC-P, ELISA. Reacts with human and mouse. For neurotransmitter biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	DBH Rabbit mAb
Remarks/Alias	DBM; ORTHYP1
Species	Human
Gene ID (Human)	1621
Gene ID	1621
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 388-617 of human DBH (NP_000778.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09172
Protein Weight	69 kDa
Shipping	Ice bag