


Anti-Dbi Rabbit Monoclonal Antibody For Wb And Elisa - P07108

Item Number: CM0005272



Introduction

Human DBI (ACBP) rabbit monoclonal antibody validated for Western blot and ELISA. Detects a 10-12 kDa band, suitable for studying acyl-CoA metabolism and GABAA receptor modulation.

[Learn More](#)

Parameter	Value
Product Name	DBI Rabbit mAb
Remarks/Alias	EP; ACBP; ACBD1; CCK-RP; DBI
Species	Human
GeneID	1622
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-104 of human DBI. (NP_001171513.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P07108
Protein Weight	10-12 kDa
Shipping	Ice bag