


Anti-Dbi Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P07108

Item Number: CM0026543



Introduction

High-affinity rabbit polyclonal antibody against human Diazepam-binding inhibitor (DBI/ACBP), validated for WB, IP, and ELISA; detects a 10-12 kDa band; cross-reacts with human DBI; suitable for lipid metabolism, neuropeptide signaling, and intracellular acyl-CoA transport research.

[Learn More](#)

Parameter	Value
Product Name	DBI Rabbit pAb
Remarks/Alias	EP; ACBP; ACBD1; CCK-RP; DBI
Species	Human
GeneID (Human)	1622
GeneID	1622
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-104 of human DBI (NP_001171513.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P07108
Protein Weight	10-12 kDa
Shipping	Ice bag