


Anti-Dbi Rabbit Polymab® For Wb, Ihc-P, Elisa - P07108 / P31786

Item Number: CM0012619



Introduction

Anti-DBI rabbit recombinant monoclonal antibody (PolymAb®) validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for acyl-CoA-binding protein/diazepam-binding inhibitor research in metabolism and neurology.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DBI Rabbit PolymAb®
Remarks/Alias	EP; ACBP; ACBD1; CCK-RP; endozepine
Species	Human/Mouse
GeneID (Human)	1622
GeneID	1622/13167
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07108/P31786
Protein Weight	10-12 kDa
Shipping	Ice bag