


Anti-Dbnl Polyclonal Antibody For Wb, Elisa - Q9Uju6

Item Number: CM0016138



Introduction

Rabbit polyclonal antibody targeting human DBNL (HIP-55, SH3P7), validated for WB and ELISA; cross-reacts with mouse and rat. DBNL functions in actin cytoskeleton reorganization, endocytosis, and T-cell activation.

[Learn More](#)

Parameter	Value
Product Name	DBNL Rabbit pAb
Remarks/Alias	ABP1; HIP55; SH3P7; HIP-55; DBNL
Species	Human
GeneID (Human)	28988
GeneID	28988
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 131-430 of human DBNL (NP_001014436.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJU6
Protein Weight	48kDa
Shipping	Ice bag