

Anti-Dbp Rabbit Polyclonal Antibody For Wb And Elisa - Q10586

Item Number: CM0017264

Introduction

High-quality rabbit polyclonal antibody against human DBP (D site-binding protein) for Western blot and ELISA applications. Detects a 34 kDa nuclear transcription factor involved in circadian rhythm regulation. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	DBP Rabbit pAb
Remarks/Alias	DABP; taxREB302; DBP
Species	Human
Gene ID (Human)	1628
Gene ID	1628
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human DBP (NP_001343.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q10586
Protein Weight	34kDa
Shipping	Ice bag