


Anti-Dbt Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip - P11182

Item Number: CM0007620



Introduction

High-quality rabbit monoclonal antibody targeting DBT (branched-chain alpha-keto acid dehydrogenase complex E2). Validated for WB, IHC-P, IP, and ELISA with human, mouse, rat cross-reactivity. Suitable for mitochondrial research and MSUD studies.

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] DBT Rabbit mAb
Remarks/Alias	E2; E2B; BCATE2; BCKADE2; BCKAD-E2; BCKDH-E2; BCOADC-E2; [KD Validated] DBT
Species	Human
Gene ID (Human)	1629
Gene ID	1629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 62-247 of DBT (NP_001909.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11182
Protein Weight	53 kDa
Shipping	Ice bag