


Anti-Dbt Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11182

Item Number: CM0022384



Introduction

Rabbit anti-DBT pAb for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Targets BCKDH E2 component; useful for branched-chain amino acid catabolism and MSUD research. P11182.

[Learn More](#)

Parameter	Value
Product Name	DBT Rabbit pAb
Remarks/Alias	E2; E2B; BCATE2; BCKADE2; BCKAD-E2; BCKDH-E2; BCOADC-E2; DBT
Species	Human
Gene ID	1629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-247 of human DBT (NP_001909.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11182
Protein Weight	53kDa
Shipping	Ice bag