


Anti-Dc-Signr/Cd299 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - Q9H2X3

Item Number: CM0006285



Introduction

Recombinant rabbit monoclonal antibody against human DC-SIGNR/CD299 (CLEC4M, L-SIGN). Validated for WB, IF/ICC, FC, ELISA. Detects a band of ~24-45 kDa. Pathogen recognition receptor involved in viral entry.

[Learn More](#)

Parameter	Value
Product Name	DC-SIGNR/CD299 Rabbit mAb
Remarks/Alias	CD299; LSIGN; CD209L; L-SIGN; DCSIGNR; HP10347; DC-SIGN2; DC-SIGNR; CD299/DC-SIGNR
Species	Human
Gene ID (Human)	10332
Gene ID	10332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 71-399 of human DC-SIGNR/CD299 (NP_055072.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	Q9H2X3
Protein Weight	24-45KDa
Shipping	Ice bag