

Anti-Dcaf12 Polyclonal Antibody For Wb And Elisa - Q5T6F0

Item Number: CM0019974



Introduction

Rabbit polyclonal antibody against human DCAF12 (DDB1- and CUL4-associated factor 12), validated for Western blot and ELISA. Recognizes human and rat DCAF12, suitable for ubiquitination and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	DCAF12 Rabbit pAb
Remarks/Alias	CT102; TCC52; WDR40A; KIAA1892; DCAF12
Species	Human
GeneID (human)	25853
GeneID	25853
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human DCAF12 (NP_056212.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q5T6F0
Protein Weight	51kDa
Shipping	Ice bag