

Anti-Dcaf7 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61962

Item Number: CM0028889

Introduction

Rabbit polyclonal antibody targeting human DCAF7 (DDB1- and CUL4-associated factor 7). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for developmental biology and ubiquitin ligase research.

[Learn More](#)

Parameter	Value
Product Name	DCAF7 Rabbit pAb
Remarks/Alias	AN11, HAN11, WDR68, SWAN-1, DCAF7
Species	Human
GeneID (Human)	10238
GeneID	10238
Immunogen	Recombinant fusion protein containing amino acids 1-342 of human DCAF7 (NP_005819.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61962
Protein Weight	39kDa
Shipping	Ice bag